





Advanced Head Rest

The newly designed advanced head rest can be purchased as an upgrade to any I-Portal® NOTC. It builds on Neuro Kinetics' unique temple stabilization system by incorporating a single mechanism that automatically centers the head and enhances stability.

This is ideal for clinical research where minimizing head position artifacts between trials is vital. Enhanced restriction is important for studies using off vertical axis rotation (OVAR) and high acceleration tests such as the NOTC Head Impulse Test.



Chin strap for added support



Auto-centering arms

- Unique auto-centering arm ensures the head is on the center of rotation.
- Temple stabilization minimizes slippage.
- Chin strap reduces slouching by providing added support.
- One button quick release allows the operator to quickly secure and release.